



Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R RPD	Blank Spike / LCS %R RPD
Sulfide, SM 4500-S2 D-0, Water (SM 4500-S2 D-00) in Water							
Preservation: ZnOAc, Cool <6°C (NaOH added at lab)							
Container: Small OJ, 500 mL, ZnOAc		Minimum Sample Volume: 250 mL		Hold Time: 7 days			
Sulfide	0.0500 mg/L	0.0500 mg/L	0.0500 mg/L		20	75 - 125	75 - 125 20