



Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R RPD	Blank Spike / LCS %R RPD
Phosphorus, Ortho-P, SM 4500-P E-1999 (SRP) (SM 4500-P E-99) in Water							
Preservation: Cool <6°C							
Container: Small OJ, 500 mL							
Minimum Sample Volume: 500 mL							
Hold Time: 2 days							
Orthophosphorus	0.00400 mg-P/L	0.00400 mg-P/L	0.00400 mg-P/L		20	75 - 125	90 - 110 20