



Analytical Method Information

Analyte	CAS #	DL	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R	Blank Spike / LCS %R	RPD
Sulfide, SM 4500-S2 D-0, Water (SM 4500-S2 D-11) in Water								
Preservation: ZnOAc, Cool <6°C (NaOH added at lab)								
Container: Small OJ, 500 mL, ZnOAc								
Minimum Sample Volume: 250 mL								
Hold Time: 7 days								
Sulfide	18496-25-8	0.0500 mg/L	0.0500 mg/L	0.0500 mg/L	20	75 - 125	75 - 125	20