



Analytical Method Information

Analyte	CAS #	DL	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R	Matrix Spike RPD	Blank Spike / LCS %R	Blank Spike / LCS RPD
Carbon, Organic Dissolved SM 5310 B (SM 5310 B-11) in Water									
Preservation: H2SO4, pH <2									
Container: Glass NM, Amber, 250 mL, Minimum Sample Volume: 250 mL Hold Time: 28 days 9N H2SO4									
Dissolved Organic Carbon		0.5000 mg/L	0.5000 mg/L	0.5000 mg/L	20	75 - 125	20	90 - 110	20
Carbon, Organic Dissolved SM 5310 B (SM 5310 B-11) in Drinking Water									
Preservation: H2SO4, pH <2									
Container: Glass NM, Amber, 250 mL, Minimum Sample Volume: 250 mL Hold Time: 28 days 9N H2SO4									
Dissolved Organic Carbon		0.5000 mg/L	0.5000 mg/L	0.5000 mg/L	20	75 - 125	20	90 - 110	20