



Analytical Method Information

Analyte	CAS #	DL	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R	Matrix Spike RPD	Blank Spike / LCS %R	Blank Spike / LCS RPD
Cyanide, Total, SM 4500-CN⁻ E-99 (SM 4500-CN⁻ E-99) in Water									
Preservation: pH>12; NaOH, Cool <6°C									
Container: Small OJ, 500 mL, NaOH									
Minimum Sample Volume: 250 mL									
Hold Time: 14 days									
Cyanide, Total	57-12-5	0.00500 mg/L	0.00500 mg/L	0.00500 mg/L	20	75 - 125	20	75 - 125	20