



### Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R RPD	Blank Spike / LCS %R RPD
<b>Cyanide, Total, SM 4500-CN<sup>-</sup> E-99 (SM 4500-CN<sup>-</sup> E-99) in Water</b>							
Preservation: pH>12; NaOH, Cool <6°C							
Container: Small OJ, 500 mL, NaOH							
Minimum Sample Volume: 250 mL							
Hold Time: 14 days							
Cyanide, Total	0.00500 mg/L	0.00500 mg/L	0.00500 mg/L		20	75 - 125	75 - 125 20