



Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R RPD	Blank Spike / LCS %R RPD
Nitrogen, Total Kjeldahl, SM 4500-Norg, Solid (SM 4500-Norg D-97) in Solid							
Preservation: Cool <6°C							
Container: Glass WM, Clear, 4 oz		Minimum Sample Weight: 100 g		Hold Time: 28 days			
Total Kjeldahl Nitrogen	41.7 mg-N/kg	41.7 mg-N/kg	41.7 mg-N/kg		20	75 - 125	90 - 110 20