



### Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R RPD	Blank Spike / LCS %R RPD
<b>Solids, Settleable SM 2540 F-97 (SM 2540 F-97) in Water</b>							
Preservation: Cool <6°C							
Container: HDPE NM, 1000 mL							
Minimum Sample Volume: 1000 mL							
Hold Time: 2 days							
Settleable Solids	0.100 mL/L	0.100 mL/L	0.100 mL/L		20		