



Analytical Method Information

| Analyte | DL | LOD | LOQ / RL | Surrogate %R | Duplicate RPD | Matrix Spike %R | RPD | Blank Spike / LCS %R | RPD |
|---------------------------------------|------------|------------|------------|-----------------|------------------|--------------------|-----|-------------------------|-----|
| TPH AK102 (AK102 DRO) in Solid | | | | | | | | | |
| Preservation: Cool <6°C | | | | | | | | | |
| Container: Glass WM, Clear, 8 oz | | | | | | | | | |
| Minimum Sample Weight: 15 g | | | | | | | | | |
| Hold Time: 14 days | | | | | | | | | |
| Diesel Range Organics (C10-C25) | 2.34 mg/kg | 2.50 mg/kg | 5.00 mg/kg | | 30 | 75 - 125 | 30 | 75 - 125 | 30 |
| surr: o-Terphenyl | | | | 50 - 150 | | | | | |