



### Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R	RPD	Blank Spike / LCS %R	RPD
<b>TPH AK102 (AK102 DRO) in Water</b>									
Preservation: Cool <6°C									
Container: Glass NM, Amber, 500 mL									
Minimum Sample Volume: 1000 mL									
Hold Time: 7 days									
Diesel Range Organics (C10-C25)	0.0340 mg/L	0.0500 mg/L	0.100 mg/L		30	75 - 125	30	75 - 125	30
surr: o-Terphenyl				50 - 150					