



Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R	RPD	Blank Spike / LCS %R	RPD
TPH AK103 (AK103 RRO) in Water									
Preservation: Cool <6°C									
Container: Glass NM, Amber, 500 mL									
Minimum Sample Volume: 1000 mL									
Hold Time: 7 days									
Motor Oil Range Organics (C25-C36)	0.0660 mg/L	0.100 mg/L	0.200 mg/L		30	60 - 120	30	60 - 120	30
surr: o-Terphenyl				50 - 150					
surr: n-Triacontane				50 - 150					