



Analytical Method Information

Analyte	DL	LOD	LOQ / RL	Surrogate %R	Duplicate RPD	Matrix Spike %R	RPD	Blank Spike / LCS %R	RPD
WA EPH Aromatics (WA EPH) in Water									
Preservation: Cool <6°C									
Container: Glass NM, Amber, 500 mL									
Minimum Sample Volume: 500 mL									
Hold Time: 7 days									
C8-C10 Aromatics	7.04 ug/L	40.0 ug/L	40.0 ug/L						
C10-C12 Aromatics	2.12 ug/L	40.0 ug/L	40.0 ug/L		30	29 - 120	30	29 - 120	30
C12-C16 Aromatics	2.55 ug/L	40.0 ug/L	40.0 ug/L		30	39 - 120	30	39 - 120	30
C16-C21 Aromatics	8.90 ug/L	40.0 ug/L	40.0 ug/L		30	54 - 120	30	54 - 120	30
C21-C34 Aromatics	14.6 ug/L	40.0 ug/L	40.0 ug/L		30	29 - 120	30	29 - 120	30
surr: o-Terphenyl					41 - 120				